**Instructions**

Bleaching gels developed solely for in-office bleaching in the dental practice. A double-cartridge for instant mixing (4:1) contains bleaching gel and activator. The mixed gel contains 6% hydrogen peroxide. In-office whitening gels may only be sold to dentists. The first round of treatment using tooth whitening products must be performed by a dentist, or under direct supervision of a dentist if equivalent safety can be guaranteed. May thereafter be handed over for continued treatment.

**Treatment periods**

PerfectWhite<sup>evo</sup> 4 x 12 min

**Indications**

- Brightening of devitalized teeth
- Brightening of pathologically discolored teeth as an alternative to crown or veneer
- Treatment of discolorations caused by pharmaceuticals (e.g. tetracycline)
- Aesthetic reasons

**Contraindications**

- Exposed dentine
- Defective filling margins
- Patients with hypersensitivities
- Allergies or incompatibilities with the ingredients of the product
- Pregnancy or breast-feeding (based upon a lack of experience)

A slight sensitization of the teeth during and after the treatment may occur in some patients. However, any discomfort normally vanishes within a relatively short period of time. **Should the patient experience excessive discomfort or pain during treatment, the bleaching process should be terminated immediately.**

**Preparation**

Determine the shade of the tooth by comparison

with a Vita<sup>®</sup> shade guide (photo 1). For assessment, re-arrange the Vita<sup>®</sup> shade guide as listed below. Values arranged in accordance with the shades: A1 A2 A3 A3.5 B1 B2 B3 B4 C1 C2 C3 C4 D1 D2 D3 D4 Values arranged in accordance with the brightness: B1 A1 B2 D2 A2 C1 C2 D3 A3 D4 B3 A3.5 B4 C3 A4 C4 A3 to B1 = increase of value by 2 A3 to A2 = increase of value by 4

Insert a cheek or lip retractor (photo 2). In case you start by cleaning the patient's teeth with air abrasion, make sure the powder is bicarbonate based. Polish teeth thoroughly with pumice without additives (photo 3). Rinse and blow completely dry. Do not let the patient rinse their mouth during treatment, as the glucoprotein in saliva significantly decreases the effect of the bleaching gel. Apply the gingiva protector along the gum line (photo 4), and cure it with the BrilliantSmile LED Plasma Light (photo 5). Do not use anesthesia. Pure activator fluid is caustic.

**The ready-mixed bleaching gel irritates the skin and the mucous membranes, and it causes severe eye damage. Always wear protective gloves and goggles when using the material! Avoid contact of the material with eyes, skin and mucous membranes. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a doctor.**

**Application**

Remove the cap from the cartridge and replace with the mixing tip. Press the gel out of the cartridge until a uniformly colored gel emerges from the mixing tip. Discard the inhomogeneous material which emerged initially. Storage at higher temperatures for a longer period of time can cause the development of isolated bubbles within the material. This may lead to inhomogeneous, liquid mixtures, especially when there is less than a quarter of the material left inside the double cartridge. In this case, do not use the remaining material. Apply the gel in a layer of 1-2 mm thickness onto the labial surfaces of the teeth to be bleached (photo 6). The mixing tip is very well suited for spreading the gel. Any gel on the gums or soft tissue must be

removed immediately. The effectiveness of the product is increased by the use of the BrilliantSmile LED Plasma Light (photo 7). **The patient and dentist should wear protective eye-wear.** Remove the gel by suction after 12 min, rinse off any residue with water and blow dry (photo 8). The process is repeated three more times. Before repeating the process, make sure the gingiva protection remains intact. If the material remaining in the cartridge is sufficient for another treatment, remove the mixing tip and close the cartridge with the blue cap. Store the material in a refrigerator. Use a new mixing tip for each patient.

**Note:** Composite fillings or ceramic and metal restorations cannot be bleached. Such fillings or restorations should be replaced about two weeks after the treatment.

**Devitalized teeth**

When treating devitalized teeth, the gel is also applied onto the lingual surface of the teeth. In order to do so, the root canal filling is removed up to about 2 mm above the enamel cement margin. Afterward, it is tightly covered with glass ionomer cement. The bleaching process is now carried out as described above.

**Concluding the treatment**

After the last bleaching process, remove the gel by suction, and thoroughly rinse off any residue with water and blow dry. Determine the tooth shade once again (photo 10). Afterward, the teeth should be remineralized with a fluoride solution or gel. Such solution or gel is also suitable for treatment of sensitization or pain occurring during the bleaching process. After treatment, the teeth are extra susceptible to discoloration. Offer the patient the BrilliantSmile NanoSeal Total+ varnish (photo 9), otherwise the patient must avoid acidic and discoloring foods, beverages and tobacco for 72 hours, in addition to rinsing with fluoride twice daily, in order to keep the best possible result.

**Shelf-life/Storage**

Gel should not be used after its expiration date, as a sufficient bleaching performance can then no longer be guaranteed. **Store the bleaching gel in a refrigerator (3-9°C), but do not freeze.** **Product contains hydrogen peroxide, potassium nitrate, sodium fluoride, sodium hydroxide, glycol derivatives, thickening agent and dye.**

**Possible sources of error**

Bleaching result insufficient	Repeat bleaching process (up to five treatment cycles in one sitting). Additional bleaching process can be performed after one week.
Patient experiences excessive discomfort or pain during treatment	Gel is in contact with exposed dentine or an untight filling (check before treatment). Patient has very sensitive teeth (depends on the individual tooth structure)
White stains on the enamel after bleaching	Local demineralization or dehydration: this is reversible and disappears after fluoridation
Gingiva appears white after bleaching	This is caused by an insufficient or untight gingiva protection, but is also reversible.
Gel protruding from the cartridge is not uniformly colored	Caused by air bubbles in the gel, which appear if the gel is stored at higher temperatures than recommended or if it is stored overlong peroxide residue
Composite enamel adhesion decreased after treatment	Wait 1-2 weeks before replacing restorations

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