i-PAK®

EN	Knitted retraction cord
LT	Gestrickte Retrakcionsschnur
DE	Gestrickte Retrakcionsschnur
BG	Ъболекарски конец
CS	Pletená navíjecí šňůra
DA	Strikket retraktionstråd
ET	Kootud retraktsiooniniit
FR	Cordon de rétraction tricoté
EL	Πλεγμένο νήμα απώθησης
HU	Kötött behúzószinór
IT	Cordone di retrazione a maglia
LV	Adīta ievilkšanas aukla
NO	Strikket retraksjonssnor
PL	Dzianinowa nić retrakcyjna
PT	Cordão de retração de malha
RO	Fir de retracție tricotat
SK	Pletená navíjacia šnúra
SL	Pletena retrakcijska nitka
ES	Cordón de retracción de punto
SV	Stickad frilågningstråd
TR	Dokuma retraksiyon kordonu
RU	Вязаная ретракционная нить
UA	В'язана ретракційна нитка

INSTRUCTION FOR USE EN

DESCRIPTION

i-PAK is non-sterile cord made of 100% cotton, knitted into thousands of tiny loops forming long interlocking chains. i-PAK is available in 4 sizes: #1, #0, #00, #000. Marking depends on the thickness of cord. A unique attribute of this cord is that it doesn't entangle in the diamond bur.

COMPOSITION

Cotton 100%. Product does not contain medicinal substance, including human blood or plasma derivative; tissues or cells, or their derivatives, of human origin; tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

Colour	knitted retraction cord 000 - black knitted retraction cord 00 - yellow knitted retraction cord 0 - lilac knitted retraction cord 1 - blue
Diameter	knitted retraction cord 000 - 1.0mm (ultra-thin) knitted retraction cord 00 - 1.2mm (thin) knitted retraction cord 0 - 1.35mm (medium) knitted retraction cord 1 - 1.5mm (medium thick)

INTENDED PURPOSE AND CLINICAL BENEFITS

i-PAK is intended for use during dental procedures for gingival retraction.

CLINICAL INDICATIONS

- **#1 for the anterior teeth preparation;**
- **#0 for anterior teeth; as upper cord for "two cord" technique;**
- **#00 during preparation and cementation of veneers; restorative procedures dealing with thin, friable tissues;**
- **#000 for lower anteriors; when luting near gingival and subgingival veneers; class III, IV V restorations; first cord for "two cord" technique.**

CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS

None known.

UNDESIRABLE SIDE EFFECTS

In susceptible individuals, product may cause allergic or irritation reactions (skin, eye, mucosa, respiratory tract irritation). Do not use product if patient is allergic to product or to any of the ingredients.

RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OR BODY FLUIDS

Part of the body – mouth. Tissues or body fluids contacted by the device – tooth, saliva, mucosa.

INTENDED USER

i-PAK is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental materials. There is no need for specific training.

STERILITY

i-PAK is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that

the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT

i-PAK IS designed to be used in dental office where ambient temperature is 18-25°C. Do not re-use.

CONSUMABLE COMPONENTS AND ACCESSORIES

No consumable components and accessories are supplied with the device.

INSTRUCTION FOR USE

1. Prepare and dry the tooth.
2. Choose the cord size. Because cord compresses upon packing, use a cord size that appears too large.
3. Cut the required length of retraction cord.
4. Soak the cord in hemostatic liquid/gel (not applicable for impregnated cord).
5. After preparing and drying the tooth, retraction cord is placed closely around the cervix and inserted carefully into the sulcus using a small retraction cord packer or spatula. The thin cord packer quickly slips it into position.
6. During the following preparation for the procedures, the mechanical effect of the cord and chemical effect of hemostatic material, results in a definite retraction of the gingiva.
7. Perform the necessary procedures.
8. If any bleeding is present, repeat steps until hemostasis is achieved.
9. When the procedures are completed, remove cord and rinse with water.

WARNINGS

Do not use product for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. Product does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS

Wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

SHelf-LIFE

Shelf-life of i-PAK is 6 years from the date of manufacture. Do not use after the expiry date. See packing for expiry date. Do not re-use.

STORAGE

Keep product tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

i-PAK is safe and performs as intended if it is used in accordance to manufacturer's instruction for use.

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

PACKAGING

REF IPU11 Knitted retraction cord #1 (244cm), medium thick, blue

REF IPU01 Knitted retraction cord #0 (244cm), medium, lilac

REF IPU02 Knitted retraction cord #00 (244cm), thin, yellow

REF IPU03 Knitted retraction cord #000 (244cm), ultra-thin, black

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