# AFFINIS<sup>®</sup> System 360

<b>DE</b>	Gebrauchsinformation	3	<b>EL</b>	Πληροφορίες χρήσης	62
<b>EN</b>	Instructions for use	8	<b>PT</b>	Instruções de Utilização	68
<b>FR</b>	Mode d'emploi	14	<b>TR</b>	Kullanma talimatı	74
<b>ES</b>	Información para el uso	20	<b>RU</b>	Инструкции по применению	80
<b>IT</b>	Informazioni per l'uso	26			
<b>SV</b>	Bruksanvisning	32			
<b>NL</b>	Gebruiksaanwijzing	38			
<b>DA</b>	Brugsanvisning	44			
<b>NO</b>	Bruksanvisning	50			
<b>FI</b>	Käyttöohje	56			

## Definition

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**Definition**

AFFINIS SYS360 is a silicone-based dental impression material indicated as a tray material for use with commercially-available automated mixing machines.

	AFFINIS® putty	AFFINIS® heavy body	AFFINIS® MonoBody
<b>Material type</b>	Polyvinylsiloxane, addition-type silicone elastomer		
<b>Base</b>	blue grey	light brown	ice blue
<b>Catalyst</b>	grey	caramel	light grey
<b>ISO</b>	ISO 4823, Type 0 putty consistency	ISO 4823, Type 1 high consistency	ISO 4823, Type 2 medium consistency
<b>Impression techniques with AFFINIS tray material</b>			
Monophase impressions			●
Wash material for the Mono-phase impressions			●
One step putty/wash impressions (double mix)	●	●	●
Two step putty/wash impressions	●	●	●
Functional impressions		●	●
Transfer coping impressions / pick-up impressions		●	●
Relining impressions		●	●
Temporary impressions	●		
Impressions for implantology (after healing is complete)	●	●	●
Copperband impressions	●		

	AFFINIS® putty	AFFINIS® heavy body	AFFINIS® MonoBody
<b>Clinical times</b>			
<b>Important</b>			
For extended working time, material can be placed in the refrigerator prior to use.	●	●	●
Direct application in the palm of the hand (additional kneading) will accelerate the setting process.	●		
Higher temperature in the surrounding environment will accelerate the setting process; lower temperatures will slow it down.	●	●	●
Check intraorally if the material is completely set before removing.	●	●	●

### Contra-indications

Provided that the product is used as directed, there are no known contra-indications.

### Side effects / Interactions

Polyvinylsiloxanes have a very good biological compatibility and up to now no harmful reactions or secondary effects on patients and/or dental personnel are known. Impression materials are intended to set in the mouth of the patient. However, they should remain in the mouth not more than twice the recommended setting time. Although they have reasonably high tear strength, care should be taken that no portion of the impression material remains in the interdental spaces or in the sulcus. Undercuts should, in certain instances, be blocked out before taking the impression.

### Trays

Tray selection depends on the impression technique and personal preference of stock trays, double arch bite trays or custom trays. For optimum adhesion, we recommend applying a thin layer of Coltène® Adhesive onto the tray, or any other brand that is indicated for use with polyvinylsiloxane impression materials.

### Preparation

1. Hold the SYS360 Cartridge facing upwards with the outlet channel closed.
2. Gently clip the safety cap of the outlet channel while holding it upwards; then remove and discard the safety cap.
3. Insert the SYS360 Cartridge according to the manufacturer's instructions for use of the mixing machine.
4. Start the mixing machine and extrude material out of the outlet

channel without the mixing tip until equal amounts of base and catalyst are dispensed. This ensures that the plungers are at the same height to achieve an optimal mix.

5. Carefully wipe off extruded material using a paper towel to avoid contamination between the base and catalyst in the outlet channels.
6. Attach a yellow dynamic mixing tip onto the outlet channel of the SYS360 Cartridge. Pay attention that the hexagon section of the yellow dynamic mixing tip is correctly aligned onto the drive shaft.
7. When the yellow dynamic mixing tip is correctly placed, the yellow fixation ring must be positioned over the dynamic mixing tip and turned right (clockwise ¼ turn) until it is securely fastened. The mixing components inside of the yellow dynamic mixing tip must rotate when the machine is started.
8. Before every use reject the first 3 cm of the mixed impression material and press out the following material without interruption.
9. Fill the mixed impression material directly into the selected impression tray. In order to prevent the formation of air voids, the yellow dynamic mixing tip should always remain immersed in the material.
10. In order to prevent excess material coming out, the plungers of the mixing machine should be disengaged after extruding material.
11. The yellow dynamic mixing tip should be left on the SYS360 Cartridge until the next application, since it acts as a seal. It should be wiped with a disinfection solution.
12. To remove the used yellow dynamic mixing tip, the fixation ring must be turned counter-clockwise and removed. The fixation ring is reusable.
13. Before the next application, check to ensure that the outlet channels of the SYS360 Cartridge are not clogged and that equal

amounts of material can be extruded. Proceed according to point 5 as follows.

14. Remove the empty SYS360 Cartridge out of the mixing machine and discard it according to local regulation.

### Notice

#### SYS360-Cartridge

- Do not drop the SYS360 Cartridge, since this may damage it, leaving it defective.

#### Yellow dynamic mixing tip

- The yellow dynamic mixing tip should be left on the SYS360 Cartridge until the next application, since it acts as a seal.
- To guarantee trouble-free application, we recommend that the Coltène/Whaledent yellow dynamic mixing tips be used exclusively.

### Recommendations for use

#### One step putty/wash technique (double mix)

While loading the tray, the dentist may begin to apply the wash-material around the preparation from the cartridge or syringe (MonoBody). Tray loading and application of wash-material around the preparation must be completed at the same time. Keep the oral tip immersed in the material at all times to avoid trapping air. *Immediately place the filled tray into the patient's mouth. Press for 2-3 s and hold in position until the material is completely set.*

#### “Two step technique”

If a 2-Step «heavy body/wash» technique is used, the primary impression made with heavy body or putty-material, must be carefully cleaned and dried before further use. This will help guarantee good adhesion between the wash- and tray-material.

### Important

Always wear gloves.

Skin secretions, latex gloves and surfaces contaminated by latex gloves can interfere with the setting process of polyvinylsiloxanes. The material and surfaces where the impression will take place (teeth, preparations, retraction cords, etc.) should only come into contact with thoroughly washed and rinsed gloves (wash 15 s with a detergent, rinse with lukewarm tap water for another 15 s) or vinyl-/nitrile gloves. Compounds containing eugenol or hemostatic substances can also impede the setting process. If (H<sub>2</sub>O<sub>2</sub>) hydrogen peroxide is used for disinfection, it is recommended to thoroughly rinse with lukewarm water to avoid the formation of air voids.

### Disinfection

The impression should be rinsed under running tap water after removal from the mouth. After rinsing, disinfection with a suitable commercial dental disinfectant solution (according to the manufacturer's recommendations) which will not affect the impression surface or dimensions. Acrylic trays should be protected against water absorption.

### Optional: autoclavable impressions

The following points should be observed when autoclaving impressions:

- Only use autoclavable components (e.g. PRESIDENT Tray AC, Coltène Adhesive AC).
- Rinse and clean impression thoroughly under flowing, lukewarm water.
- Impressions can be autoclaved directly after taking the impression.
- Autoclave in steam sterilizer at 134°C/ 273°F (Prion program)

When sterilizing implant impressions, be sure clarify in advance with the manufacturer, whether or not the implant components (e.g. impression posts, etc.) are autoclavable.

#### **Fabrication of models**

The stone models can be poured after 30 min at the earliest. The impression remains dimensionally stable for a practically unlimited period of time (for at least 7 days). The surface tension can be reduced and pouring can be facilitated when the impression is briefly washed out with a detergent and rinsed thoroughly in lukewarm clear water afterwards. All commercially-available industry-standard dental stone model materials (i.e. Fujirock Dental Stone, Hard Rock Dental Stone) can be used.

#### **Electro plating**

The impression is compatible with all commercially available copper- and silver plating baths.

#### **Cleaning of trays**

The impression can be removed with a blunt instrument. The Coltène® Adhesive can be removed by soaking it in a universal, commercially-available solvent or petroleum ether. Solvents should only be used in a well-ventilated area. Clean and disinfect trays according to customary methods. PRESIDENT disposable impression trays are designed for single-use and are therefore not reusable.

#### **Shelf life and storage**

AFFINIS SYS360 will fulfil its intended purpose at least until the expiry date, which is displayed on the containers, when kept in well-sealed containers at temperatures 15-23 °C / 59-73 °F and 50 % relative air humidity. The SYS360 Cartridge should be stored in a horizontal posi-

tion. Impressions should be stored at normal room temperature. Avoid exposure to heat and direct exposure to sunlight.

#### **Caution**

Federal law restricts this device to sale by or on the order of a dentist.

#### **Marking**

The expiry date and [LOT]-number are displayed on the package and cartridge.

#### **Technical data ISO 4823:2000**

Measurements are based on 23 °C / 73 °F room temperature and 50% relative humidity.

#### **putty**

Mixing time (15ml):	0:15 min
Time in mouth (35 °C / 95 °F):	3:30 min

#### **heavy body**

Mixing time (15ml):	0:15 min
Total working time:	0:45 min
Time in mouth (35 °C / 95 °F):	2:00 min

#### **MonoBody**

Mixing time (15ml):	0:15 min
Total working time:	0:45 min
Time in mouth (35 °C / 95 °F):	2:30 min

#### **Date of issue**

07-2013

Only supplied to dentists and dental laboratories or upon their instructions.

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Coltène/Whaledent AG