



8280-N
to
8294-N

MEDICOM® DRY-BACK DENTAL BIBS 3-PLY

FEATURES & BENEFITS


- The Medicom® Dry-back 3-ply dental bibs range is made of 2 layers of embossed cellulose and 1 layer of waterproof polyethylene material.
- This allows maximum protection with horizontal embossed pattern and superior absorption thanks to additional layer of cellulose.



Recommended use :

- This bib, placed polyethylene face on the patient, permits to establish a barrier against the risks of infection during the dental care.

Available in 13 colours:



PRODUCT DESCRIPTION

Product name: Dry-Back Dental Bibs 3-Ply
 Product type: Single use, non-sterile
 Material: Embossed cellulose, Polyethylene
 Colour : Lavender, Blue, White, Yellow, Green, Dusty Rose, Peach, Aqua, Dark blue, Orange, Black, Lime, Turquoise
 Grammage : 33 g/m² (± 10%)
 Folding: M
 Size: 33 x 45,5 cm (± 0,5 cm)
 Unit weight: 6 g (± 10%)
 Packaging: Bulk of 500 units
 Origin: France



INDUSTRIES

Medical / Dental

8280-N
to
8294-N

MEDICOM® DRY-BACK DENTAL BIBS 3-PLY

CERTIFICATION

Manufactured under an ISO 9001 certified quality system.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C

Shelf life of the product: 5 years

LOGISTICAL INFORMATION



Outer case specifications



EAN Codes

| Item | Size mm | Gross weight kg | QTY/ carton | QTY/ pallet | Item | Outer case |
|--------|-----------------|-----------------|--------------------|-------------------------------------|--------|-------------------|
| 8280-N | | | | | 8280-N | 8 718 007 006 230 |
| 8282-N | | | | | 8282-N | 8 718 007 006 377 |
| 8283-N | | | | | 8283-N | 8 718 007 006 490 |
| 8284-N | | | | | 8284-N | 8 718 007 006 605 |
| 8285-N | | | | | 8285-N | 8 718 007 006 711 |
| 8286-N | | | | | 8286-N | 8 718 007 006 810 |
| 8288-N | 300 x 240 x 340 | 3,2 | Bulk of 500 pieces | 78 cartons (6 layers of 13 cartons) | 8288-N | 8 718 007 006 933 |
| 8289-N | | | | | 8289-N | 8 718 007 007 039 |
| 8290-N | | | | | 8290-N | 8 718 007 007 138 |
| 8291-N | | | | | 8291-N | 8 718 007 007 183 |
| 8292-N | | | | | 8292-N | 8 718 007 012 422 |
| 8293-N | | | | | 8293-N | 8 718 007 013 719 |
| 8294-N | | | | | 8294-N | 8 718 007 906 141 |



Medicom

Groupe KOLMI HOPEN SAS

Bld de la chanterie | 49124 Saint Barthélemy d'Anjou – France

Bld de la chanterie | BP 10059 – 49181 Saint Barthélemy d'Anjou Cedex – France

T : + 33 (0) 241 963 434 | F : + 33 (0) 241 963 453 | sales@kolmi-hopen.com | www.medicom-eu.com