

## TECHNICAL DATA ENBIO®

Class	B
Power supply	220 - 230V AC/ 50/60 Hz
Installed power	2250W
Max. power consumption	15A
Working pressure	2,1Bar
Max pressure	2,3 Bar
Max process temperature	137°C
Process chamber volume	2,7l
Tray inner diameter	28 x 18 x 2,8 cm
Programs	134°C FAST, 134°C, 121°C
Max. packages amount (circa)	8
Max. weight of the load	0,5 kg
Mass	15 kg
Level of security	IP20
Data archiving	USB Memory stick
Water supply*	external

\*demineralized or distilled water only

